

Monthly radioactivity measurement "March, 2012" was conducted.**Measurement Result:**

No radioactive substances (Iodine/Caesium) was detected with any of Products/Cannulae.

Table 1. Radioactivity measurement - March/2012

Sample Category	Nuclide					
	Iodine	Caesium			Radiation Dose	
	I-131	Cs-134	Cs-137	Cs-136		
Unit	Bq/kg(L)	Bq/kg(L)	Bq/kg(L)	Bq/kg(L)	µSv/h	cpm
Products (Needles)	ND ^{*1)}	ND ^{*1)}	ND ^{*1)}	ND ^{*1)}	-	-
Cleaning Water	Measurement once in two months (Next April)					
Cannulae	-	-	-	-	ND ^{*2)}	-
Components	Measurement once in two months (Next April)					

*1) Not Detected (Detection Limit=7)

*2) Not Detected (Equivalent count rate as Background)