## Monthly radioactivity measurement "April, 2012" was conducted.

## **Measurement Result:**

No radioactive substances (Iodine/Caecium) was detected with any of Products/Cannulae.

Table 1. Radioactivity measurement - April/2012

, , , , , , , , , , , , , , , , , , , ,						
Sample Category	Nuclide					
	Iodine	Caesium			Dadiation Dage	
	I-131	Cs-134	Cs-137	Cs-136	Radiation Dose	
Unit	Bq/kg(L)	Bq/kg(L)	Bq/kg(L)	Bq/kg(L)	μSv/h	cpm
Products (Needle)	ND <sup>*1</sup>	ND <sup>*1</sup>	ND <sup>*1</sup>	ND <sup>*1</sup>	-	-
Cleaning Water	ND <sup>*1</sup>	ND <sup>*1</sup>	ND <sup>*1</sup>	ND <sup>*1</sup>	-	-
Cannulae	Measurement once in two months (Next May)					
Components		_	-	_	ND <sup>*2</sup>	-

<sup>\*1)</sup> Not Detected (Detection Limit=9)

<sup>\*2)</sup> Not Detected (Equivalent count rate as Background)